

QUICK VIEW

Cleaning System Information



F3900

Do's

- Follow the manufacturer cleaning instructions on the cleaning powder bottle.
- Gently brush your dental appliance after the 15 minute cleaning cycle to rinse and remove any loosened debris.
- **Use (1) ultrasonic cleaning cycle during the (15-minute) appliance soaking period as a daily cleaning regimen.** For occasional deep cleaning if needed; use (2) - (15-minute) appliance soakings with (2) ultrasonic cleaning cycles (rinse in-between soaking periods).
- Always pour used cleaning solution out from the front of the ultrasonic unit/tank, and pour slowly to avoid splashing which can get into your eyes or ruin clothing.

Dont's

- Do not soak dental appliances longer than the recommended 15 minutes.
- Do not use cleaning powder on dental appliances containing soldered joints, as they will become tarnished.
- Deep cleaning is not necessary if you are following the recommended daily cleaning regimen.
- Do not pour used cleaning solution out from the side or back of the ultrasonic unit. Water will damage the electronics.
- Do not run more than (2) ultrasonic cycles consecutively or you will overheat and damage the transducer.
- Do not use the plastic tray that comes with the unit when cleaning, because it reduces the ultrasonic cleaning power.

CAUTION: Please remember that you are using a professional strength cleaning powder that is not to be compared to over-the-counter dental appliance cleaning tablets. Keep away from children.

Please refer to your Cleaning Powder Handling Instructions.

For any questions regarding your ultrasonic cleaning unit call: iSonic at 1 (847)-850-0404